English INSTRUCTION FOR USE (IFU)

Sterile, Single Use
UNICAPRONE
Monofilament Polyglecaprone 25
Synthetic Absorbable Surgical
SUTURES / LIGATURES

DESCRIPTION: UNICAPRONE (Polyglecaprone 25) is a sterile, synthetic, absorbable, monofilament surgical suture manufactured from Poly (glycolide-co-ε-caprolactone). UNICAPRONE sutures are either un-dyed (appears Beige in colour) or dyed with D&C violet # 2 (Colour index number 60725) during polymerization. UNICAPRONE sutures are available in numerous gauge sizes and lengths, non-needed or attached to stainless steel needles. UNICAPRONE sutures have uniform diameter with clear surface, low stiffness and have the property of smooth passage through tissue. Polyglecaprone 25 sutures are relatively inert, non-antigenic, non-pyrogenic and elicit only a mild tissue reaction during absorption. UNICAPRONE complies with the requirements of the European Pharmacopoeia for sterile synthetic absorbable monofilament sutures.

INDICATIONS: UNICAPRONE sutures are indicated for use in general soft tissue approximation and/or ligation where an absorbable material is indicated but not for use in cardiovascular or neurological tissues, in microsurgery or ophthalmic surgery.

APPLICATION: Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

PERFORMANCE: UNICAPRONE (Polyglecaprone 25) elicits a minimal initial inflammatory reaction in tissue. Progressive loss of tensile strength and eventual absorption of UNICAPRONE occurs by means of hydrolysis where the polymer degrades to adipic acid which is subsequently absorbed and metabolized in the body. Absorption begins as loss of tensile strength followed by loss of mass. Absorption studies show the following profile for degradation of the suture. All of the original tensile strength is lost by 21 days and the absorption is essentially complete in 120 days.

Post Implantation	Approx. Original Strength Remaining
7 days - undyed suture	50-60%
7 days - dyed suture	50-60%
14 days - undyed suture	20-30%
14 days - dyed suture	25-35%

CONTRAINDICATIONS: These sutures, being absorbable should not be used where extended approximation of tissues under stress is required.

WARNINGS/PRECAUTIONS/INTERACTIONS: Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing UNICAPRONE suture for wound closure, as risk of wound dehiscence may vary with the site of application and the type of suture used. Physicians should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. As an absorbable suture, UNICAPRONE suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention, or which may require additional support. Skin sutures which remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. UNICAPRONE suture knots must be properly placed to be secure. Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be

particularly appropriate when knotting monofilaments. Avoid prolonged exposure to elevated temperatures. Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping at the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in "sharps" containers.

ADVERSE REACTIONS: Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

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