English INSTRUCTION FOR USE (IFU)

Sterile, Single Use
UNICRYL quick
Braided Coated Polyglycolic Acid (PGA)
(Fast Absorbing) Synthetic Absorbable Surgical SUTURES

DESCRIPTION: UNICRYL quick is a synthetic quicky absorbable surgical suture material composed of a homopolymer of glycolide. The braid is coated with a mixture of polycaprolactone and calcium stearate. UNICRYL quick sutures are supplied un-dyed (appears Beige in colour). The characteristic of quick loss of strength is achieved by use of a polymer material with a lower molecular weight than normal polyglycolic acid suture.

Although this suture is a synthetic absorbable suture, its performance characteristics are intended to model the performance of collagen (surgical gut) suture. The knot tensile strength of UNICRYL quick suture meets EP knot tensile strength requirements for collagen sutures.

UNICRYL quick is available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes. UNICRYL quick complies with the requirements of the European Pharmacopoeia for sterile braided synthetic absorbable surgical sutures and the USP for synthetic absorbable sutures.

INDICATIONS: UNICRYL quick suture is indicated only for use in superficial soft tissue approximation of the skin and mucosa, where only short term wound support (7-10 days) is required. UNICRYL quick suture is not intended for use in ligation, ophthalmic, cardiovascular or neurological procedures.

APPLICATION: Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

PERFORMANCE: UNICRYL quick suture, when used in closure of skin and mucous membranes, typically begins to fall off 7-10 days post-operatively and can be wiped off subsequently with sterile gauze. Natural mechanical abrasion of the sutures while in situ may also accelerate this disappearance rate. Rapid loss of tensile strength may preclude the need for sitch removal. UNICRYL quick elicits a minimal to moderate acute inflammatory reaction in tissue.

Absorption begins as a loss of tensile strength followed by a loss of mass. Progressive loss of tensile strength and eventual absorption occurs by means of hydrolysis, where the homopolymer degrades to glycolic acid, which is subsequently absorbed, and metabolized by the body.

Absorption studies of UNICRYL quick sutures show that approximately 50% of the original tensile strength remains during 5-8 days. All of the original tensile strength is lost by approximately 14 to 18 days and the absorption is approximately complete by 42 days.

CONTRAINDICATIONS: Due to the rapid loss of tensile strength, this suture should not be used where extended approximation of tissues under stress is required or where wound support beyond 7 days is required.

WARNINGS / PRECAUTIONS / INTERACTIONS: Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing UNICRYL quick suture for wound closure, as risk of wound dehiscence may vary with the type of application and the type of suture used. Physicians should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions, which may delay wound healing. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, UNICRYL quick suture may act transiently as a foreign body.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention, or which may require additional support. Skin sutures which remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as

suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. UNICRYL quick suture, which is treated with coating to enhance handling characteristics, requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. Avoid prolonged exposure to elevated temperatures. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

ADVERSE REACTIONS: Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distention occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local contact with the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

STERILITY (STERILIZATION): UNICRYL quick sutures are sterilized by ethylene oxide gas. Single use only. Do not re-sterilize. Do not use if package is opened or damaged. Discard open unused sutures.

STORAGE: Recommended storage conditions: below 25°C, away from moisture and direct heat. Do not use after the expiry date.

INDICATIONS: UNICRYL quick suture is indicated only for use in superficial soft tissue approximation of the skin and mucosa, where only short term wound support (7-10 days) is required. UNICRYL quick suture is not intended for use in ligation, ophthalmic, cardiovascular or neurological procedures.

APPLICATION: Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

PERFORMANCE: UNICRYL quick suture, when used in closure of skin and mucous membranes, typically begins to fall off 7-10 days post-operatively and can be wiped off subsequently with sterile gauze. Natural mechanical abrasion of the sutures while in situ may also accelerate this disappearance rate. Rapid loss of tensile strength may preclude the need for sitch removal. UNICRYL quick elicits a minimal to moderate acute inflammatory reaction in tissue.

Absorption begins as a loss of tensile strength followed by a loss of mass. Progressive loss of tensile strength and eventual absorption occurs by means of hydrolysis, where the homopolymer degrades to glycolic acid, which is subsequently absorbed, and metabolized by the body.

Absorption studies of UNICRYL quick sutures show that approximately 50% of the original tensile strength remains during 5-8 days. All of the original tensile strength is lost by approximately 14 to 18 days and the absorption is approximately complete by 42 days.

CONTRAINDICATIONS: Due to the rapid loss of tensile strength, this suture should not be used where extended approximation of tissues under stress is required or where wound support beyond 7 days is required.

WARNINGS / PRECAUTIONS / INTERACTIONS: Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing UNICRYL quick suture for wound closure, as risk of wound dehiscence may vary with the type of application and the type of suture used. Physicians should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions, which may delay wound healing. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, UNICRYL quick suture may act transiently as a foreign body.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention, or which may require additional support. Skin sutures which remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as

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