



SURGICAL GLOVES

ORTHO POWDER FREE

SURGICAL GLOVES

- POWDER FREE
- MICRO ROUGHENED
- NON SLIP SURFACE FOR OPTIMAL SENSITIVITY
- THICK MATERIAL WITH HIGH TENSILE STRENGTH
- LOW PROTEIN CONTENT
- NATURAL RUBBER LATEX
- NATURAL COMFORT
- STERILE
- MANUFACTURED IN ACCORDANCE WITH EN 455

Product Code

SOPF55
SOPF60
SOPF65
SOPF70
SOPF75
SOPF80
SOPF85
SOPF90

Product Description

Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 5.5
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 6.0
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 6.5
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 7.0
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 7.5
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 8.0
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 8.5
Sterile Single Use Natural Rubber Latex Orthopaedic Surgical Gloves Powder Free 9.0

UNIMED[®]
GLOVES

Factory of Medical Rubber and Plastic Materials
(A Branch of United Medical Industries Co. Ltd).

www.unimed.com.sa

CAT/SOPF Rev.1

Issue Date: - 01/07/2017