



UNICRYL

3 Weeks
Effective Tensile
Strength

UNIMED[®]
SUTURES


United Medical Industries Co. Ltd.

www.unimed.com.sa

UNICRYL

- Coated Braided Polyglycolic Acid B.P. Violet and Undyed.
- Synthetic absorbable sterile surgical suture composed of a homopolymer of glycolide.
- Available in range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- Minimal tissue reaction.
- Easy passage through tissue.
- Secure Knot tying and easy knot tie down.
- Intended for use where an absorbable suture or ligature is indicated. It can be used in general soft tissue approximation or ligation and in ophthalmic procedures.
- Sterilised by Ethylene Oxide.
- Coating – Polyglycolic acid is a homopolymer of glycolic acid, which is braided into sutures of various sizes. The braid is coated with mixture of Calcium Stearate and Polycaprolactone or Poly-D-Lactic acid and N-e-lauroyl-L-lysine.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United State Pharmacopoeia (U.S.P.).

UNICRYL - Approximate Breaking Strength Retention



Head Office

P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel. +966 1 2650665
Fax. +966 1 2650569
general@unimed.com.sa